

2 Start new Sensor with Reader

STEP 1



Press Home Button to turn on Reader.


Note: If using Reader for the first time, follow the prompts to set date and time.

STEP 2



Touch **Start New Sensor**.
Hold Reader within 1.5 inches (4 cm) of Sensor to scan it.

STEP 3



Review the important information on the screen.

Sensor can be used to check your glucose after 60 minutes.

3 Check your glucose

STEP 1




Press Home Button to turn on Reader.

STEP 2



Hold Reader within 1.5 inches (4 cm) of Sensor to scan it.

STEP 3



Reader shows your glucose reading. This includes your Current Glucose, Glucose Trend Arrow, and Glucose Graph.

Setting Alarms

The Sensor automatically communicates with the Reader and can give you glucose alarms. Alarms are on by default. To change their settings or turn alarms off, follow these steps.

IMPORTANT: Glucose alarms are an important safety feature. Please consult your health care professional before making changes.

STEP 1



Press Home Button to go to Home Screen. Touch .

STEP 2



Touch **Alarms** and then touch **Change Alarm Settings**.

STEP 3




Select and set your alarms. Touch **done** to save.

Using Alarms

You must scan your Sensor to see your glucose reading.



Touch **Dismiss Alarm & Check Glucose** or press the Home Button to dismiss alarm and check your glucose.



Touch **No** to dismiss alarm.
Touch **Yes** or press the Home Button to dismiss alarm and scan the Sensor.



When any alarm is on, sound/vibration and signal status display on the Home Screen.

| Symbol | What it means |
|--------|--------------------------------------|
| | Sound and Vibration ON |
| | Sound and Vibration OFF |
| | Sensor communicating with Reader |
| | Sensor not communicating with Reader |

Refer to User's Manual for additional symbols.

If you have followed the instructions described in the User's Manual and are still having difficulty setting up your System or if you are not sure about a message or reading, contact your health care professional.

The shape of the circle sensor unit, FreeStyle, Libre, and related brand marks are owned by Abbott.

©2019 Abbott. 487463792-001 Rev. C (9/20)

Consult instructions for use

Manufacturer
Abbott Diabetes Care Inc.
1360 South Loop Road
Alameda, CA 94502 USA

Abbott